



MAK-1900 CITY WATER/BOTTLE DELIVERY SYSTEM WITH AIR PURGE

INSTRUCTIONS FOR USE

*Please ensure you have read and followed the installation instructions prior to operation.

PRODUCT DESCRIPTION

Treat your waterlines while remaining on city water...with the flip of a switch.

The MAK-1900 bottle system is exclusively designed for dental offices that use city water and want to treat and maintain clean dental unit waterlines. Using an ingeniously simple method of dual toggle switches, the MAK-1900 allows for regular waterline cleaning in any operator that delivers city water directly into the dental unit without a bottle system.

FOR GENERAL DENTAL PROCEDURES

- 1). After complete installation, flip Toggle switch marked CITY WATER to ON.

MAK-1900 CITY WATER/BOTTLE DELIVERY SYSTEM WITH AIR PURGE

INSTRUCTIONS FOR USE

DIRECTIONS FOR WATERLINE CLEANING

1. Turn dental unit OFF. Remove bottle and fill with Mint-A-Kleen. Place bottle back onto dental unit. Turn dental unit ON.
2. Turn CITY WATER toggle switch to OFF.
3. Turn DUWL (Dental Unit Waterline) CLEANING toggle switch to AIR.
4. Air purge waterlines by activating handpiece hoses and air/water syringes until water is eliminated.*
5. Turn DUWL CLEANING toggle to MINTAKLEEN.
6. Activate handpiece hoses and air/water syringes until Mint-A-Kleen exits.
7. Turn dental unit OFF. Allow Mint-A-Kleen to treat waterlines for at least 12 hours (preferably over the weekend before proceeding to next step).
8. Turn dental unit ON. Turn CITY WATER TOGGLE SWITCH to ON. Flush waterlines with water for 30 seconds to remove any excess Mint-A-Kleen.
9. Continue with dental procedures.

* Note – For optimal cleaning, air purge and clean all instruments and tubing that contain water using Mint-A-Kleen and its protocol. Additionally, clean any areas that cannot be reached with normal protocol which are called 'Dead-Ends'. These dead-ends can be opened and cleaned with a Drain Tube. A Drain Tube is included with the MAK-1900.

* Activate all the water outlets including the quick-disconnect water outlet for the cavitron until air comes out. Use the Drain Tube to open quick-disconnect water outlet for the cavitron. Connect Drain Tube to the saliva ejector first and then to the quick-disconnect water outlet. After air purging is completed, disconnect the Drain Tube.